



**All Advanced Infusion (AITECS) Syringe Pumps include the following features:**

- Ideal for all anesthesia infusion applications, including **propofol** and **fentanyl**
- Exclusive standby feature
- Parameter retention
- High Flow rates 1500 mL/h
- Bolus protection feature
- Integral I V Pole Clamp
- One Year Warranty — Extended Warranty available
- In-Stock ready for immediate shipment
- Volume Discounts



Wilburn Medical  
USA

1-877-WILBURN  
1-877-945-2876

Wilburn Medical USA  
146 Furlong Industrial Dr.  
Kernersville, NC 27284  
WilburnMedicalUSA.com

# 12S Pro

# \$2,295



- Anesthesia Syringe Pump
- Manual Bolus Mode
- Automatic Bolus Mode
- Automatic Dose Rate Calculation
- µg/kg/min Parameter
- Incremental & Decremental Programming
- Preset Drug Profiles

# 12S Auto

**With Profile Mode**

# \$2,695



- Anesthesia Syringe Pump
- Manual Bolus Mode
- Automatic Bolus Mode
- Automatic Dose Rate Calculation
- µg/kg/min Parameter
- Incremental & Decremental Programming
- Preset Drug Profiles

FEATURE	BENEFIT
Unique syringe confirmation system.	Ensures proper flow rate for each infusion regardless of change in syringe manufacturers.
Unique operator programmable occlusion selection system.	Eliminates needless nuisance alarms caused by positional IV lines.
Exclusive operator programmable standby feature.	Pre-program your pump to prepare for your next patient and you can now put the pump on hold. Eliminate noisy and needless standby alarm completely.

**WilburnMedicalUSA.com** **1-877-WILBURN**  
**(1-877-945-2876)**

WilburnMedicalUSA.com is an extensive online medical supply catalog offering tens of thousands of items as well as detailed information about our most popular product lines. This site also features monthly specials, discounts, and clearance items. Shop online with your credit card or use our innovative Quote List to request quotes online.

## SyringeInfusion.com

SyringeInfusion.com is an extension of WilburnMedicalUSA.com which is dedicated to the many lines of syringe pumps we offer, with a particular focus on the Aitecs line of syringe pumps. Visit SyringeInfusion.com and view our two fully interactive Flash tutorials for the Aitecs syringe pumps.